

# Cadd-Legacy<sup>®</sup> 1 Ambulatory Infusion Pump

## CADD LEGACY PUMPS ARE DESIGNED TO PROVIDE:

- CADD<sup>®</sup> pump reliability
- Single-therapy pump simplicity
- Enhanced programming flexibility
- Ease of use

Infusion rates from 1 to 3,000 ml/day can meet the requirements of a wide range of continuous delivery protocols.

## DESIGNED FOR CONTINUOUS INFUSION THERAPIES

Delivery Modes	Continuous
Delivery Routes	Intravenous, Intra-arterial, Subcutaneous, Intraperitoneal, Epidural, Intrathecal
Programming Units	ml
Reservoir Volume	1 ml–9,999 ml or not in use
Continuous Rate	1–3,000 ml/24 hr

### Easy To Read Display

- Provides information for programming adjustments and troubleshooting

### Cassette detection

- Designed to detect cassette attachment

### Occlusion detection

- Upstream and downstream occlusion sensors alert the patient and clinic

### Air-in-line detection

- Off, High- and Low-sensitivity settings is designed to alert the patient and clinician if air is present in the administration set tubing



### Axelerist Inc.

80 Business Park Dr., Ste. 209  
Armonk, NY 10504  
United States  
(914) 253-8840  
[info@axelerist.com](mailto:info@axelerist.com)

### Axelerist B.V.

Westplein 12 12-14  
3016BM Rotterdam  
The Netherlands  
0031 (0)10 30 70 217  
[info@axelerist.com](mailto:info@axelerist.com)