

BD Alaris™ 8015 PC

The BD Alaris™ 8015 PC point-of-care unit provides a common user interface for programming the 8100 AND 8110 infusion modules, which aids in reducing the complexity at the point of care. The BD Alaris™ 8015 PC is intended for adult, pediatric, and neonatal care in hospital and health care facilities.



Dimensions	6.9" W x 8.8" H x 9" D
Weight	7.2 lbs
Battery Operation	Depends on the number of modules attached and module activity - Approximately 6 hours with 1 pump module - Approximately 3 hours with 4 pump modules - Wireless connectivity and internal backup lithium battery for data storage
Communication Data Port	RS-232 with an RJ45 connector
Power Requirements	100-240V, 50/60 Hz, 150VA Max

Biomed Tested / Calibrated

Axelerist Inc.

80 Business Park Dr., Ste. 209
Armonk, NY 10504
United States
(914) 253-8840
info@axelerist.com

Axelerist B.V.

Westplein 12 12-14
3016BM Rotterdam
The Netherlands
0031 (0)10 30 70 217
info@axelerist.com



Equipping You
for Success.

www.axelerist.com