

BD Alaris™ 8100

The BD Alaris™ 8100 infusion pump is a large volume pump (LVP) that can help protect every infusion type, including continuous and intermittent deliveries. It is designed for the delivery of fluids, medications, and blood products for adult, pediatric, and neonatal patients.



Dimensions	8" H x 5.5" D x 3.3" W
Weight	2.5 lbs
KVO Rate	1 ml/hour
Flow Rates	0.1 -99.9 mL/hour in 0.1 mL increments 100-999 mL/hour in 1 mL increments
Maximum Infusion Pressure	654 mmHg

Biomed Tested / Calibrated

Axelerist Inc.

80 Business Park Dr., Ste. 209
Armonk, NY 10504
United States
(914) 253-8840
info@axelerist.com

Axelerist B.V.

Westplein 12 12-14
3016BM Rotterdam
The Netherlands
0031 (0)10 30 70 217
info@axelerist.com



Equipping You
for Success.

www.axelerist.com