

# BD Alaris™ 8110

The BD Alaris™ 8110 is a syringe module that combines precision instrumentation with advanced software to help ensure the right dose rate. It is intended for the delivery of fluids, medications, blood, and blood products using intermittent or continuous delivery. This syringe pump module uses a pressure sensing disc which shortens the time to alarm and reduces the delays upon infusion start up.



|                  |   |
|------------------|---|
| Dimensions       | 4.5" W x 15" H x 7.5" D   |
| Weight           | 4.5 lbs   |
| Flow Rates       | 0.01-9.99 ml/hr in 0.01 ml increments<br>10-99.9 ml/hr in 0.1 ml increments<br>100-999 ml/hr in 1 ml increments |
| Maximum Pressure | 16 psi without pressure sensing disc<br>21 psi with pressure sensing disc                                       |
| Alerts           | Syringe unlocked<br>Infusion near end<br>High pressure/occlusion<br>Plunger disengaged<br>Low battery           |

*Biomed Tested / Calibrated*

#### **Axelerist Inc.**

80 Business Park Dr., Ste. 209  
Armonk, NY 10504  
United States  
(914) 253-8840  
[info@axelerist.com](mailto:info@axelerist.com)

#### **Axelerist B.V.**

Westplein 12 12-14  
3016BM Rotterdam  
The Netherlands  
0031 (0)10 30 70 217  
[info@axelerist.com](mailto:info@axelerist.com)



Equipping You  
for Success.

[www.axelerist.com](http://www.axelerist.com)