



# Countertop Freezer Free-Standing

## (48L) -20 °C

### BUILT FOR SCIENTIFIC AND MEDICAL APPLICATIONS

Units feature a digital microprocessor temperature controller allowing for the precise temperature management necessary for critical samples and supplies. A forced air directional refrigeration system provides superior temperature uniformity and recovery after door openings.

|                      |                            |
|----------------------|----------------------------|
| Max Capacity:        | 48 L ( 1.7 cu. ft.)        |
| Temperature:         | -15° to -25°C              |
| Alarm:               | High/Low temperature alarm |
| Exterior Dimensions: | 23.75W x 24D x 20.375H"    |
| Defrost:             | Manual                     |
| Door Lock:           | Key Lock Upper Center      |
| Voltage:             | 115V                       |
| Frequency:           | 60Hz                       |



#### **Axelerist Financial**

80 Business Park Dr., Ste. 209  
Armonk, NY 10504  
United States  
(914) 253-8840  
[info@axelerist.com](mailto:info@axelerist.com)

#### **Axelerist B.V.**

Westplein 12 12-14  
3016BM Rotterdam  
The Netherlands  
0031 (0)10 30 70 217  
[info@axelerist.com](mailto:info@axelerist.com)