

EBA 200 Ambient Centrifuge

Compact performance with built-in versatility.

The EBA 200 is a compact clinical centrifuge offering performance and safety for small sample volumes. The EBA 200 features a built-in, versatile, 8-place, fixed-angle rotor and adjustable LCD display for convenient and quick operation. Available in low and high speed (up to 8,000 RPM) models.

Adjustable parameters for flexible application control.

The EBA 200 and EBA 200 S feature a flexible control panel allowing you to adjust your application parameters. Adjust time, speed (RPM), RCF and braking with the easy to use LCD display.

Visual and audible end-of-cycle notification



Max Capacity	8 x 15 ml
Max. RPM / RCF	6,000 / 3,461
Max. Blood Tubes*	8 (up to 10 ml)
Max. Conical Tubes*	4 (15 ml)
Temperature Control	Ambient
Dims. (H x W x D)	9 x 10 x 14 in.

EBA 200 Blood Tube Package-8



Tube Qty.	Volume (mL)	Size (mm)	200		200 S	
			RPM	RCF	RPM	RCF
8	1.6-5	13x75	≤6,000	≤8,000	≤2,697	≤4,794
8	4-7	13x100	≤6,000	≤8,000	≤3,461	≤6,153
8	8.5-10	16x100	≤6,000	≤8,000	≤3,461	≤6,153
4	8	16x125	≤6,000	≤8,000	≤3,461	≤6,153

Includes: (1) Centrifuge, (1) Built-in Angle Rotor, (8) Adapters 13x75mm, (8) Adapters 13x100mm

Axelerist Inc.

80 Business Park Dr., Ste. 209
Armonk, NY 10504
United States
(914) 253-8840
info@axelerist.com

Axelerist B.V.

Westplein 12 12-14
3016BM Rotterdam
The Netherlands
0031 (0)10 30 70 217
info@axelerist.com