



# Eppendorf 5702R Refrigerated Centrifuge with A-4-38 Rotor

## COMPACT SPACE-SAVING DESIGN AND LOW NOISE LEVEL

The 5702 family of low speed, general purpose centrifuges was especially designed to meet the requirement of low- to mid-throughput clinical research and cell culture labs.

They are compact in design so that one can fit comfortably on every lab bench, quiet in operation to benefit your work environment.

Speed:	3,000 x g (4,400 rpm)
Max Capacity:	12 x 15 ml
Temperature settings:	from -9 to +40°C
Rotor Type:	Swing Out
Power Supply:	120 V, 60Hz and 230V 50/60Hz
Certifications:	230V model is CE



### Axelerist Financial

80 Business Park Dr., Ste. 209  
Armonk, NY 10504  
United States  
(914) 253-8840  
[info@axelerist.com](mailto:info@axelerist.com)

### Axelerist B.V.

Westplein 12 12-14  
3016BM Rotterdam  
The Netherlands  
0031 (0)10 30 70 217  
[info@axelerist.com](mailto:info@axelerist.com)