

High Performance Undercounter Refrigerator

+2 to +8°C Undercounter Refrigerator | 3.0 cu. ft. | 85L

This high performance undercounter refrigerator is designed to meet current and evolving CDC, NSF, and ANSI guidelines for vaccine storage. This unit is ENERGY STAR® Certified and capable of maintaining the 5°C setpoint standard for storing vaccines with a +/- uniformity of 3°C to better reduce the risk of freezing vital vaccines and biologics. Built with accessibility in mind, this unit is ADA compliant and engineered with a glass door for easy visibility of stored contents.



Exterior Dimensions	22.5" W x 22.1" D x 29.6" H 572 x 562 x 751 mm)
Interior Dimensions	18.5" W x 14.6" D x 18.9" H (471 x 370 x 480 mm)
Effective Capacity	3.0 cu.ft. (85 L)
Shelves And Shelf Dimensions	Vinyl coated wire, 2 shelves (adjustable) Dimensions: 17.3" x 14.0" (440 x 355 mm)
Access Port	1: Bottom of back of product (0.6" (16 mm) diameter)
Temperature Control Range	2 to 8 in 0.1° increments
Alarms	Power failure (B - Momentary), high temperature deviation (V-B), low temperature deviation (V-B), door open (V-B), diagnostics (V) V=Visual, B=Buzzer
Net Weight	114 lbs (52 kg)
Electrical Requirements	115V, 60Hz, 1Ø, NEMA 5-15P requires NEMA 5-15R receptacle
Doors	1 glass, self closing
Temperature Display	LED
Door Lock	Key lock

Axelerist Inc.

80 Business Park Dr., Ste. 209
Armonk, NY 10504
United States
(914) 253-8840
info@axelerist.com

Axelerist B.V.

Westplein 12 12-14
3016BM Rotterdam
The Netherlands
0031 (0)10 30 70 217
info@axelerist.com



Equipping You
for Success.

www.axelerist.com