


MIR Spirolab Spirometer

All-in-one Spirometer with 7" display, Embedded Printer and Oximetry option, to carry on the go. Desktop, Stand-alone and PC-based Spirometer, with Oximetry Option.



| | |
|---------------------|---|
| Flow sensor | bi-directional digital turbine |
| Volume rate | 10 L |
| Flow range | ±16L/s |
| Volume accuracy | ±2.5% or 50 mL |
| Flow accuracy | ±5% or 200 mL/s |
| Dynamic resistance | <0.5 cm H ₂ O/L/S |
| Temperature sensor | semiconductor (0-45°C) |
| Test available | FVC, VC, IVC, MVV, PRE-POST |
| Measured parameters | FVC, FEV1, FEV1/FVC%, FEV1/PEF, FEV1/VC, FEV1/FEF0.5, DTPEF |
| Memory capacity | Up to 10000 tests |

Axelerist Inc.

80 Business Park Dr., Ste. 209
Armonk, NY 10504
United States
(914) 253-8840
info@axelerist.com

Axelerist B.V.

Westplein 12 12-14
3016BM Rotterdam
The Netherlands
0031 (0)10 30 70 217
info@axelerist.com