



# Memory Traceable Monitoring Thermometer

Temperature Range:	-50 to 70°C (-58 to 158°F)
Temperature Resolution:	1°
Temperature Accuracy:	±1°C
Cable Length:	10'
Temp MIN/MAX Memory:	Y
Data Logging:	N
Exportable Data (Y/N):	N
LCD Digit Size:	1/2 in
Dimensions:	2-3/4" x 3/4" x 4-1/4"
CE Marked:	Yes



## Alarm:

Alarm is programmable in 1° increments, sounds continuously when temperature rises above/falls below set points and can be turned off manually

## Bottle Probe:

Bottle insulates sensor from rapid temperature changes when refrigerator door is opened. Ten-foot micro-cable permits refrigerator doors to close on it.

## Traceable® Certificates:

Multi-point calibration on an individually-numbered Traceable® Certificate which assures accuracy from our ISO/IEC 17025:2005 (1750.01) calibration laboratory accredited by A2LA. It indicates traceability of measurements to the SI units through NIST or other recognized national measurement institutes (NMI) that are signatories to the CIPM Mutual Recognition Agreement.

## Accurate and reliable:

Range is -58 to 158°F and -50 to 70°C with a resolution of 1° and accuracy of ±1°C. Ten-foot micro-cable permits refrigerator doors to close on it.

### Axelerist Financial

80 Business Park Dr., Ste. 209  
Armonk, NY 10504  
United States  
(914) 253-8840  
[info@axelerist.com](mailto:info@axelerist.com)

### Axelerist B.V.

Westplein 12 12-14  
3016BM Rotterdam  
The Netherlands  
0031 (0)10 30 70 217  
[info@axelerist.com](mailto:info@axelerist.com)