



# Mortara Eli 150rx EKG Machine

## COMPACT, LIGHTWEIGHT MACHINE WITH LEADING CLINICAL RESEARCH FUNCTIONALITY

The ELI 150 Rx is a compact and lightweight 12-lead EKG machine with leading clinical research functionality. It is capable of capturing, viewing, transmitting, printing, and storing ECG data. Internal storage allows for up to 150 digital ECG records to be saved. An easy-to-use alphanumeric elastomer keypad includes dedicated "one-touch" buttons.



Paper Type:	Z-fold thermal paper; 108 mm (4") wide, perforated double Z-fold. 200 double Z-fold sheets.
Thermal Recorder:	Computer-controlled dot array; 8 dots/mm
Thermal Writer Speeds:	5, 10, 25, or 50 mm/s
Gain Settings:	5, 10, or 20 mm/mV
ECG Storage:	Up to 150 ECGs
Weight:	7.2 lb (3.3 kg) including battery (without paper)
Dimensions:	11.25" x 11.5" x 3.75" (29.2 cm x 30.5 cm x 10.2 cm)
Power Requirements:	Universal AC power supply (100 - 240 VAC at 50/60 Hz) 50 VA; internally rechargeable battery.

### Axelerist Financial

80 Business Park Dr., Ste. 209  
Armonk, NY 10504  
United States  
(914) 253-8840  
[info@axelerist.com](mailto:info@axelerist.com)

### Axelerist B.V.

Westplein 12 12-14  
3016BM Rotterdam  
The Netherlands  
0031 (0)10 30 70 217  
[info@axelerist.com](mailto:info@axelerist.com)