

# NDD EasyOne Air Portable and PC Spirometer

EasyOne Air is a flexible, accurate and proven spirometry solution. Easily connect to the PC with Bluetooth for real-time patient incentives, data exchange and EMR connectivity.



Test Types	FVC, FVL, SVC, MVV and Pre/Post testing
Languages:	English, Danish, Dutch, French, German, Italian, Polish, Portuguese, Russian, Spanish, Swedish
EMR connectivity:	EasyOne Connect
Dimensions:	3.4" x 6.1" x 1.4", 13 oz 87 mm X 155 mm x 36 mm, 0.356 kg
Standards:	ATS/ERS 2019 and 2005, NIOSH, OSHA, SSA Disability
Features	<ul style="list-style-type: none"> <li>• Large color screen</li> <li>• Easy navigation</li> <li>• Calibration-free technology</li> <li>• Rechargeable battery: &gt;100 tests one cycle</li> <li>• Bluetooth connectivity</li> <li>• Flexible use</li> <li>• Portable or PC mode</li> <li>• Real-time animated incentives</li> <li>• Real-time data transfer</li> <li>• Proven EMR integration with top EMR/EHR systems</li> </ul>

## Axelerist Inc.

80 Business Park Dr., Ste. 209  
Armonk, NY 10504  
United States  
(914) 253-8840  
[info@axelerist.com](mailto:info@axelerist.com)

## Axelerist B.V.

Westplein 12 12-14  
3016BM Rotterdam  
The Netherlands  
0031 (0)10 30 70 217  
[info@axelerist.com](mailto:info@axelerist.com)