

# B. Braun Perfusor Space Syringe Pump

The compact, light-weight Perfusor® Space® syringe pump provides an easy solution to complex therapies with an intuitive and user-friendly design. It is an accurate, stackable syringe pump, perfect for any care setting. Available in both wireless and non-wireless options.

## Simple User Interface

- Same clean design on every pump
- Similar workflow on all Space Pumps reduces training required
- Safety prompts and reminders
- Easy to backlit screen

## Compact Design

- Small, light-weight
- Designed for easy transport within facility

## Keyguard, Non-Numeric Keypad

- Eliminates zero vs. decimal errors
- Eliminates double key bounce errors



Type of Unit	Syringe Driven Infusion Pump
Dimensions	249 x 68 x 152 mm (WxHxD), drive parked
Weight	Approx. 1.4 kg
Moisture protection	IP 22, drip protected for horizontal usage
Display	Backlit graphic display, ~40° read angle from all sides
Keypad	Backlit keys, cell phone like cursor navigation
Flow Rates	0.01-1800 ml/h
Accuracy of set delivery rate	<< ± 0,5 % mechanical accuracy ± 2 % in compliance with IEC/EN 60601- 2-24
Operating Temperature	+5° C ... + 40° C, +41° F .... +105 °F
Voltage	11 - 16 V DC supplied by external Space Power Supply or by SpaceStation
Battery operating time	Minimal power consumption of the devices and new battery technology ensure long operating times, e.g. 16.8 hrs @ 1 ml/h, 16 hrs @ 5 ml/h, 14 hrs @ 20 ml/h

## Axelerist Inc.

80 Business Park Dr., Ste. 209  
Armonk, NY 10504  
United States  
(914) 253-8840  
[info@axelerist.com](mailto:info@axelerist.com)

## Axelerist B.V.

Westplein 12 12-14  
3016BM Rotterdam  
The Netherlands  
0031 (0)10 30 70 217  
[info@axelerist.com](mailto:info@axelerist.com)