



Personal Low Temperature Incubator (18 L)

SAVE SPACE WITH A PERSONAL LOW TEMPERATURE INCUBATOR IDEAL FOR SMALL SAMPLE WORK.

Volume :	18 L
Temperature Range (°C):	17°C-40°C
Ext. Dimensions:	10.2W x 16.3H x 18.5D "
Temperature Uniformity:	± 1.2°C
Packaging:	24W x 24H x 32D "
Int. Dimensions:	7.1W x 12.2H x 11.4D "
Electrical:	100-240V AC, 50/60Hz
Temperature Setting Accuracy:	± 0.2°C
Certifications:	CE



- Monitor samples through transparent door with internal light
- Easy-to-use control panel
- Digital display for easy reading
- Common applications include microbiology, blood sample storage, bacterial research, histology, coliform determination, yeast growth.

Axelerist Financial

80 Business Park Dr., Ste. 209
Armonk, NY 10504
United States
(914) 253-8840
info@axelerist.com

Axelerist B.V.

Westplein 12 12-14
3016BM Rotterdam
The Netherlands
0031 (0)10 30 70 217
info@axelerist.com