



Portable Ultra-Low Temperature Freezer

IDEAL FOR CLINICAL TRIALS

Storage Volume:	25 liters (0.9 cu.ft.)
Temperature:	-86°C to -20°C @ 32°C (90°F) ambient
Alarm:	Warm and Cold Alarms Setpoint $\pm 10^{\circ}\text{C}$
Exterior Dimensions:	18.6" W x 19.9" D x 25.0" H (472 x 505 x 635 mm)
Net Weight, Empty:	21 kg (46 lbs.), nominal
Voltage Range, Universal:	120V or 240V, AC; 60Hz or 50Hz; or 12VDC from mobile source
Certification/Listing:	CE Certified
Storage Capacity:	18 standard 2" boxes
Security:	Lockable Lid



- The benefits of Stirling Ultracold freezers reach all levels of the clinical trial.
- Lightest overall weight, lowest shipping costs, easy to handle
- Can plug into any outlet, eliminating the cost of converting an outlet to 240V
- Smallest footprint in its class, giving the investigator more options for placement within the site.

Axelerist Financial

80 Business Park Dr., Ste. 209
Armonk, NY 10504
United States
(914) 253-8840
info@axelerist.com

Axelerist B.V.

Westplein 12 12-14
3016BM Rotterdam
The Netherlands
0031 (0)10 30 70 217
info@axelerist.com