



# SECA 217 Stadiometer

Quick and easy installation without wall mounting. The seca 217 mobile stadiometer is characterized by a particularly high quality of material and a sophisticated plug system. This gives it a uniquely good stability. The adjustable wall spacer also ensures a secure hold without any screw connection. The stadiometer is particularly suitable for mobile use in conventional medical examinations and measuring patients at home. After all, it is easy to transport the stadiometer when disassembled into its individual parts.



Net Weight	7.9 lbs (3.6 kg)
Product Width	12.9 inch (328 mm)
Product Height	84.4 inch (2,145 mm)
Product Depth	22.6 inch (574 mm)
Measuring Range	8 inch - 6 ft 9 inch (20 - 205 cm)

**Axelerist Inc.**

80 Business Park Dr., Ste. 209  
Armonk, NY 10504  
United States  
(914) 253-8840  
[info@axelerist.com](mailto:info@axelerist.com)

**Axelerist B.V.**

Westplein 12 12-14  
3016BM Rotterdam  
The Netherlands  
0031 (0)10 30 70 217  
[info@axelerist.com](mailto:info@axelerist.com)