

# SECA 869 Floor Scale with Remote Display (550LB Capacity)

Cable remote display for easy reading and operation.

The cable remote display can be mounted on the wall or set up on the table. This is how the seca 869 adapts perfectly to your workplace. You can switch on the TARE, HOLD and BMI function conveniently using the buttons next to the display. With its handle and high stability, the seca 869 flat scales are also ideal for mobile use.



Capacity	550 lbs (250 kg)
Graduation	0.2 lb < 330 lbs > 0.5 lb (100 g < 150 kg > 200 g)
Net Weight	9.3 lbs (4.2 kg)
Product Width	12.9 inch (328 mm)
Product Height	84.4 inch (2,145 mm)
Product Depth	22.6 inch (574 mm)

#### **Axelerist Inc.**

80 Business Park Dr., Ste. 209  
Armonk, NY 10504  
United States  
(914) 253-8840  
[info@axelerist.com](mailto:info@axelerist.com)

#### **Axelerist B.V.**

Westplein 12 12-14  
3016BM Rotterdam  
The Netherlands  
0031 (0)10 30 70 217  
[info@axelerist.com](mailto:info@axelerist.com)