

# Thermo Scientific Medifuge Small Benchtop Centrifuge

The Thermo Fisher Medifuge centrifuge with its unique Dualspin rotor allows for both fixed angle or swing bucket applications.

Includes DualSpin rotor body, fixed angle buckets (set of 8), swinging buckets (set of 8), green/short spacers (set of 8), yellow/long spacers (set of 8).

Three selectable deceleration profiles, including brake-off option, are designed to optimize separation.

Certifications/Compliance	UL, CE, RoHS, CSA, IVD, FDA class 1
Drive System	Induction drive
Height (English) Exterior Lid Open	20.04 in.
Length (English) Exterior	17.7 in.
Temperature Range	2°C to 40°C
Width (English) Exterior	12.8 in.
Capacity	8 x 15mL (Fixed Angle), 8 x 10 mL Swinging Bucket)
Max. Speed	4,900 rpm
Display	Large, Bright LED
Includes	DualSpin rotor body, fixed angle buckets (set of 8), swinging buckets (set of 8), green/short spacers (set of 8), yellow/long spacers (set of 8)
Profile (Acceleration/Braking)	1 accel (standard)/3 decel (standard, soft, brake-off)
Safety Features	Lid lock, steel armored ring, view port for speed calibration, lid emergency release
Max. RCF	3,114 x g (fixed angle); 3,490 x g (swinging bucket)



#### **Axelerist Inc.**

80 Business Park Dr., Ste. 209  
Armonk, NY 10504  
United States  
(914) 253-8840  
[info@axelerist.com](mailto:info@axelerist.com)

#### **Axelerist B.V.**

Westplein 12 12-14  
3016BM Rotterdam  
The Netherlands  
0031 (0)10 30 70 217  
[info@axelerist.com](mailto:info@axelerist.com)