

Traceable Big-Digit Memory Thermometer, 1 Bullet Probe

- Jumbo digital display is easy to read from distances
- Includes NIST-Traceable certificate from our A2LA accredited (1750.01) calibration laboratory
- F/°C selectable
- High/Low alarm
- Features a bullet probe



Meter Measurement Type	Digital Indicator
Number Of Probe Channels	1
Min Temperature (° F)	-58
Max Temperature (° F)	158
Min Temperature (° C)	-50
Max Temperature (° C)	70
Accuracy	±2.0°F (±1.0°C)
Resolution	0.1°
Alarm	Yes
Mounting Type	Desktop or wall-mount
Display Type	1.25" Display, 3.5-digit LCD
Height (in)	4 1/2
No. of Probe(s) Included	1
Probe Diameter (in)	0.125
Cable Length (in)	120
Battery	AAA Battery (1)

Axelerist Inc.

80 Business Park Dr., Ste. 209
Armonk, NY 10504
United States
(914) 253-8840
info@axelerist.com

Axelerist B.V.

Westplein 12 12-14
3016BM Rotterdam
The Netherlands
0031 (0)10 30 70 217
info@axelerist.com